



Use Word Finder

Use Word Finder to find a word and see where it is.

1. Select **Menu → Word Finder**.
2. Type the word, then select **Find**.
3. One or more paths will appear if the word is available.
4. Select a path to go to the word. Follow the highlighted buttons to step through the button sequence.

Edit a Button

1. Find the button you want to edit, then select **Menu → Edit Page**.
2. Select the button to edit → **Edit this Button**.
3. Change the **Button Label** to update what you see on the button.
4. Change the **Button Message** to update what the button says.
5. Change how the words on this button are pronounced:
 - a. Type how the words sound in **Alternate Pronunciation**.
 - b. Select **Pronounce** to test the message.
6. Change the **Image** for the button.
 - a. Choose **Select Image** to look for a symbol.

- b. Choose **Import Photo** to use a photo from the gallery.
- c. Choose **Use Camera** to take a photo.
- d. Choose **Clear Image** to delete the image.

7. Select **Save → Done**.

Menu Lock

Use Menu Lock to block access to the menu. There are two ways to lock the menu:

Use Menu Passcode

1. Select **Menu → Settings**. Scroll to the **Editing** section.
2. Type a passcode and toggle on **Use Menu Passcode**.
3. Select the **blue checkmark** ✓.
4. To unlock the menu, select **Menu → enter the passcode**.

If you forget the menu passcode, enter **72584556**

Remove Editing from the Navigation Bar

1. Leave the **Unity AAC** app, go to **iOS Settings**.
2. Scroll on the left panel and select **Apps**.
3. Scroll on the right panel and select **Unity AAC**.
4. Toggle **Allow Editing** to **Off**.

5. Go back to the **Unity AAC** app.
6. To Allow Editing again, go to **iOS Settings → Apps → Toggle Allow Editing On**.
7. Go back to the TouchChat app.

Change the Voice

1. Select **Menu → Settings**. Scroll to the **Speech** section.
2. Select the language/voice combination you want to change.
3. Select **Voice** to open the list of voices in that language.
4. Select the **play icon** next to a voice to try it.
 - a. Select the **voice name** to choose it if already downloaded.
 - b. Select the **download icon** then select the voice name to choose it. Wi-Fi is needed to download the voice to your iOS device.
5. Select the **back arrow** → select the **blue checkmark** ✓.

Change the Way the Device Pronounces a Word

1. Select **Menu → Edit Pronunciation**.
2. Select the **plus sign** +.
3. Type the correct spelling under **Word**.



4. Type how the word sounds under **Pronounce As**.
5. Select **Pronounce** to test it.
6. Select **Save** → the **blue checkmark** .

Vocabulary Builder

Use Vocabulary Builder to show fewer words on a page for a short time.

1. Select **Menu** → **Vocab Builder**.
2. Turn the **On/Off** toggle to ON.
3. Type the word or phrase you want to add in the “Add new word” field.
4. Select a word to add to the list.
5. Continue adding words to the list.
6. Select **Done**.
7. When you are done, select **Menu** → **Vocab Builder**.
8. Turn the **On/Off** toggle to OFF.

Managing Vocabulary Builder Lists

Save a List

1. Select **Menu** → **Vocabulary Builder** → **Save** → Collapse the keyboard. .
2. Select “Exported Text” to rename → **Save** → **Done**.

3. If you add more words to the list, make sure to save again.

Load a List

1. Select **Menu** → **Vocabulary Builder** → **Load**.
2. Choose the drive/folder your list is saved in.

Connecting Multiple Language Vocabulary Files

1. Find the button you want to use to go to another vocabulary. Select **Menu** → **Edit Page**.
2. Select the button to edit → **Edit this Button/Create New Button**.
3. Change the **Button Label**, **Button Message**, and **Image**.
4. Scroll to **Button Actions** → **Go to Vocabulary**.
5. Choose the vocabulary file you want to go to → **Save** → **Done**.
6. Select the button to go to the other vocabulary file.
7. Follow the same steps to create a button to go back to the first vocabulary file.

Back Up a Vocabulary File

Choose the backup option that works best for you: AirDrop®, Dropbox®, Google Drive®, or Email. For more specific information on using your preferred storage space, use the Unity AAC App Manual or Video Tutorials.

1. Select **Vocab** → **Choose New Vocab**.
2. Select **Menu** → **Import/Export Vocab**.
3. Select **Export and Share**.
4. Select the vocabulary file you want to back up/export.
5. Select **Done** and then select the method you want to use to share the file.

Find Video Tutorials

[PRC-Salttillo YouTube Channel](#) → **Unity AAC App** section

