

# PRC-Salttillo Accessibility Conformance Report

## Revised Section 508 Edition

### Via Pro 11

(Based on VPAT® Version 2.5)

**Name of Product/Version:** Via Pro 11

**Report Date:** June 23, 2025

**Product Description:** Via Pro is an electronic speech generating device that augments communication for an individual with speech/language impairment. The device and its language programs give the user the power to initiate conversation, seek information, state opinions, and share feelings. The display has a an 11-inch screen. Via devices are built on an iOS platform and offer the choice of three communication applications: LAMP Words for Life®, TouchChat®, or Dialogue® AAC.

**Contact Information:** Phone: (800) 262-1984 / E-mail Address: [info@prc-salttillo.com](mailto:info@prc-salttillo.com)

**Web site Address:** <https://www.prc-salttillo.com/co-branded-products>

**Notes:** Device Version: Via Pro 11.2

**Evaluation Methods Used:** Manual Review, Testing is based on general product knowledge

# Applicable Standards/Guidelines

This report covers the degree of conformance for the following accessibility standard/guidelines:

Standard/Guideline	Included In Report
<a href="#">Web Content Accessibility Guidelines 2.0</a>	Level A (Yes / No ) Level AA (Yes / No ) Level AAA (Yes / No )
<a href="#">Web Content Accessibility Guidelines 2.1</a>	Level A (Yes / No ) Level AA (Yes / No ) Level AAA (Yes / No )
<a href="#">Web Content Accessibility Guidelines 2.2</a>	Level A (Yes / No ) Level AA (Yes / No ) Level AAA (Yes / No )
<a href="#">Revised Section 508 standards published January 18, 2017 and corrected January 22, 2018</a>	(Yes / No )
<a href="#">EN 301 549 Accessibility requirements for ICT products and services - V3.1.1 (2019-11)</a> AND <a href="#">EN 301 549 Accessibility requirements for ICT products and services - V3.2.1 (2021-03)</a>	(Yes / No )

## Terms

The terms used in the Conformance Level information are defined as follows:

- **Supports:** The functionality of the product has at least one method that meets the criterion without known defects or meets with equivalent facilitation.
- **Partially Supports:** Some functionality of the product does not meet the criterion.
- **Does Not Support:** The majority of product functionality does not meet the criterion.
- **Not Applicable:** The criterion is not relevant to the product.
- **Not Evaluated:** The product has not been evaluated against the criterion. This can only be used in WCAG Level AAA criteria.

# WCAG 2.0 Report

Tables 1 and 2 also document conformance with Revised Section 508:

- Chapter 5 – 501.1 Scope, 504.2 Content Creation or Editing
- Chapter 6 – 602.3 Electronic Support Documentation

Note: When reporting on conformance with the WCAG 2.0 Success Criteria, they are scoped for full pages, complete processes, and accessibility-supported ways of using technology as documented in the [WCAG 2.0 Conformance Requirements](#).

## Table 1: Success Criteria, Level A

Notes:

Criteria	Conformance Level	Remarks and Explanations
<b><u>1.1.1 Non-text Content</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Each representative button on the device has a graphic image that is a representation of contextual text for that button. Synthesized prompts are included for text alternatives for users to identify graphic images. Authoring Tool: Not Applicable
<b><u>1.2.1 Audio-only and Video-only (Prerecorded)</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: The Device provides all audio with representative text. Authoring Tool: Not Applicable
<b><u>1.2.2 Captions (Prerecorded)</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Software: Not Applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Device software does not contain any prerecorded audio or video. Authoring Tool: Not Applicable
<b><u>1.2.3 Audio Description or Media Alternative (Prerecorded)</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Not applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Device does not have video or produce time-based media. Authoring Tool: Not Applicable
<b><u>1.3.1 Info and Relationships</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports	Web: Not Applicable Electronic Docs: Not Applicable

Criteria	Conformance Level	Remarks and Explanations
<ul style="list-style-type: none"> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Authoring Tool: Not Applicable	Software: Information, structure, and relationships conveyed through presentation can be programmatically determined or are available in text. Authoring Tool: Not Applicable
<b><u>1.3.2 Meaningful Sequence</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Device Software supports sequential information through word power features. Authoring Tool: Not Applicable
<b><u>1.3.3 Sensory Characteristics</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Operating device content does not rely solely on sensory characteristics of components such as shape, size, visual location, orientation, or sound. Authoring Tool: Not Applicable
<b><u>1.4.1 Use of Color</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Color used on graphic images and buttons are not the only visual means of conveying information, indicating an action, prompting a response, or distinguishing a visual element. Authoring Tool: Not Applicable
<b><u>1.4.2 Audio Control</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Device supports audio volume control. Authoring Tool: Not Applicable
<b><u>2.1.1 Keyboard</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports	Web: Not Applicable Electronic Docs: Not Applicable Software: A keyboard is built within the device for the user to utilize.

Criteria	Conformance Level	Remarks and Explanations
<ul style="list-style-type: none"> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Authoring Tool: Not Applicable	Authoring Tool: Not Applicable
<b><u>2.1.2 No Keyboard Trap</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Software: The keyboard is built within the software does not allow traps. Authoring Tool: Not Applicable
<b><u>2.2.1 Timing Adjustable</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Device supports user adjustable timing. Authoring Tool: Not Applicable
<b><u>2.2.2 Pause, Stop, Hide</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Device supports user adjustable controls. Authoring Tool: Not Applicable
<b><u>2.3.1 Three Flashes or Below Threshold</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Device supports Three flashes or below threshold. Authoring Tool: Not Applicable
<b><u>2.4.1 Bypass Blocks</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software) – Does not apply to non-web software</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs) – Does not apply to non-web docs</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Authoring Tool: Not Applicable

Criteria	Conformance Level	Remarks and Explanations
<b><u>2.4.2 Page Titled</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable
<b><u>2.4.3 Focus Order</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable
<b><u>2.4.4 Link Purpose (In Context)</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Device supports that the purpose of each link can be determined from the link text alone or from the link text together with its programmatically determined link context. Authoring Tool: Not Applicable
<b><u>3.1.1 Language of Page</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable
<b><u>3.2.1 On Focus</u></b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Software: Device Supports when any component receives focus, it does not initiate a change of context. Authoring Tool: Not Applicable
<b><u>3.2.2 On Input</u></b> (Level A) Also applies to:	Web: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable

Criteria	Conformance Level	Remarks and Explanations
Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Software: Changing settings does not change context. Authoring Tool: Not Applicable
<b>3.3.1 Error Identification</b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Electronic Docs: Software: Supports Authoring Tool:	Web: Electronic Docs: Software: Descriptive texts are provided to explain through input errors. Authoring Tool:
<b>3.3.2 Labels or Instructions</b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Labels and instructions are provided for device functions with the appropriate manual. Authoring Tool: Not Applicable
<b>4.1.1 Parsing</b> (Level A) Also applies to: WCAG 2.0 – Always answer ‘Supports’ Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Application is not software code requiring parsing. Authoring Tool: Not Applicable
<b>4.1.2 Name, Role, Value</b> (Level A) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: : The devices use standard controls for accessible technologies. Authoring Tool: Not Applicable



## Table 2: Success Criteria, Level AA

Notes:

Criteria	Conformance Level	Remarks and Explanations
<a href="#">1.2.4 Captions (Live)</a> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable
<a href="#">1.2.5 Audio Description (Prerecorded)</a> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable
<a href="#">1.4.3 Contrast (Minimum)</a> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: The device conforms to the visual presentation of text and images of text have an adequate contrast ratio. Authoring Tool: Not Applicable
<a href="#">1.4.4 Resize text</a> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Device supports resize of text. Authoring Tool: Not Applicable
<a href="#">1.4.5 Images of Text</a> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports	Web: Not Applicable Electronic Docs: Not Applicable Software: Device supports images of text. Authoring Tool: Not Applicable

Criteria	Conformance Level	Remarks and Explanations
<ul style="list-style-type: none"> <li>602.3 (Support Docs)</li> </ul>	Authoring Tool: Not Applicable	
<b>2.4.5 Multiple Ways</b> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software) – Does not apply to non-web software</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs) – Does not apply to non-web docs</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Authoring Tool: Not Applicable	Web: Device is not a website. Electronic Docs: Not Applicable Authoring Tool: Not Applicable
<b>2.4.6 Headings and Labels</b> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Device supports proper labels for user elements. Authoring Tool: Not Applicable
<b>2.4.7 Focus Visible</b> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Keyboard indicator is focus visible. Authoring Tool: Not Applicable
<b>3.1.2 Language of Parts</b> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Supports Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Device supports language of parts. Authoring Tool: Not Applicable
<b>3.2.3 Consistent Navigation</b> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software) – Does not apply to non-web software</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs) – Does not apply to non-web docs</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Authoring Tool: Not Applicable

Criteria	Conformance Level	Remarks and Explanations
<a href="#">3.2.4 Consistent Identification</a> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software) – Does not apply to non-web software</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs) – Does not apply to non-web docs</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Authoring Tool: Not Applicable
<a href="#">3.3.3 Error Suggestion</a> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable
<a href="#">3.3.4 Error Prevention (Legal, Financial, Data)</a> (Level AA) Also applies to: Revised Section 508 <ul style="list-style-type: none"> <li>501 (Web)(Software)</li> <li>504.2 (Authoring Tool)</li> <li>602.3 (Support Docs)</li> </ul>	Web: Not Applicable Electronic Docs: Not Applicable Software: Not Applicable Authoring Tool: Not Applicable	Web: Not Applicable Electronic Docs: Not Applicable Software: Not applicable, the device does not utilize legal, financial or data information. Authoring Tool: Not Applicable

### Table 3: Success Criteria, Level AAA

Notes: **Not Evaluated:** The product has not been evaluated against the criterion.

Revised Section 508 – Level AAA Does not apply

Criteria	Conformance Level	Remarks and Explanations
<a href="#">1.2.6 Sign Language (Prerecorded)</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">1.2.7 Extended Audio Description (Prerecorded)</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">1.2.8 Media Alternative (Prerecorded)</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">1.2.9 Audio-only (Live)</a> (Level AAA)	Web:	Web:

Criteria	Conformance Level	Remarks and Explanations
Revised Section 508 – Does not apply		
<a href="#">1.4.6 Contrast (Enhanced)</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">1.4.7 Low or No Background Audio</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">1.4.8 Visual Presentation</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">1.4.9 Images of Text (No Exception)</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">2.1.3 Keyboard (No Exception)</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">2.2.3 No Timing</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">2.2.4 Interruptions</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">2.2.5 Re-authenticating</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">2.3.2 Three Flashes</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">2.4.8 Location</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">2.4.9 Link Purpose (Link Only)</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">2.4.10 Section Headings</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">3.1.3 Unusual Words</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">3.1.4 Abbreviations</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">3.1.5 Reading Level</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">3.1.6 Pronunciation</a> (Level AAA)	Web:	Web:

Criteria	Conformance Level	Remarks and Explanations
Revised Section 508 – Does not apply		
<a href="#">3.2.5 Change on Request</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">3.3.5 Help</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:
<a href="#">3.3.6 Error Prevention (All)</a> (Level AAA) Revised Section 508 – Does not apply	Web:	Web:

# Revised Section 508 Report

Notes:

## Chapter 3: [Functional Performance Criteria \(FPC\)](#)

Notes:

Criteria	Conformance Level	Remarks and Explanations
302.1 Without Vision	Partially Supports	Without other supports, vision is required to operate the touchscreen on the device. Auditory Prompting and Auditory Fishing are provided for auditory cue when scanning or using direct selection with the touchscreen. The device supports users with low vision; the device is not designed for users without any vision.
302.2 With Limited Vision	Partially Supports	Some functionality is supported for limited vision by enlarging font size, text color and background colors in the main text box (via the device settings), to better accommodate high contrast viewing and color blindness for persons with limited vision.
302.3 Without Perception of Color	Supports	Device supports operation without Perception of color.
302.4 Without Hearing	Supports	Device supports operation without hearing.
302.5 With Limited Hearing	Supports	Device supports operation with limited hearing.
302.6 Without Speech	Supports	Device supports operation without speech.
302.7 With Limited Manipulation	Supports	Device supports operation with limited manipulation by on-screen elements and menus on the Device. Also, by way of one

Criteria	Conformance Level	Remarks and Explanations
		or more optional attachment switch controller(s).
302.8 With Limited Reach and Strength	Supports	Device supports operation with limited reach and strength.
302.9 With Limited Language, Cognitive, and Learning Abilities	Supports	Device supports operation with limited language, cognitive, and learning abilities.

## Chapter 4: [Hardware](#)

Notes:

Criteria	Conformance Level	Remarks and Explanations
<a href="#">402 Closed Functionality</a>	Heading cell – no response required	Heading cell – no response required
<b>402.1 General</b>	Heading cell – no response required	Heading cell – no response required
<b>402.2 Speech-Output Enabled</b>	Heading cell – no response required	Heading cell – no response required
402.2.1 Information Displayed On-Screen	Supports	Information is clearly displayed on-screen.
402.2.2 Transactional Outputs	Supports	Although there is no transactional Outputs associated with this device by default as expressed by the VPAT guidelines, this criterion is considered supported.
402.2.3 Speech Delivery Type and Coordination	Supports	Device supports Speech Delivery Type and Coordination.
402.2.4 User Control	Supports	Device supports user control via touchscreen to activate on-screen elements and menus.
402.2.5 Braille Instructions	Does Not Support	This device does not provide braille instructions.
<b>402.3 Volume</b>	Heading cell – no response required	Heading cell – no response required
402.3.1 Private Listening	Supports	Device Supports Private Listening.
402.3.2 Non-private Listening	Supports	Device supports Non-private Listening.
402.4 Characters on Display Screens	Supports	Device supports characters on the display screen as it relates to volume control.

Criteria	Conformance Level	Remarks and Explanations
402.5 Characters on Variable Message Signs	Not Applicable	Device is not an Aisle Seats, Designated Aisle Seats. Aisles, Parking, Passenger Loading zones, Assembly areas, Checkout aisles.
<b><u>403 Biometrics</u></b>	Heading cell – no response required	Heading cell – no response required
403.1 General	Not Applicable	Device does not utilize biometrics.
<b><u>404 Preservation of Information Provided for Accessibility</u></b>	Heading cell – no response required	Heading cell – no response required
404.1 General	Supports	Device supports preservation of Information transmitted.
<b><u>405 Privacy</u></b>	Heading cell – no response required	Heading cell – no response required
405.1 General	Supports	Device supports the same degree of privacy of input and output.
<b><u>406 Standard Connections</u></b>	Heading cell – no response required	Heading cell – no response required
406.1 General	Supports	Device supports at least one of each type of industry standard non-proprietary format connections.
<b><u>407 Operable Parts</u></b>	Heading cell – no response required	Heading cell – no response required
407.2 Contrast	Supports	Device supports contrast control.
<b><u>407.3 Input Controls</u></b>	Heading cell – no response required	Heading cell – no response required
407.3.1 Tactilely Discernible	Supports	Device supports tactile discernibility with “Keyguard” or “TouchGuide”. Various settings are available to prevent accidental activations including an option to activate upon release (as opposed to activating on touch) as well as dwell and acceptance time setting to accommodate those that can’t come off the keys as easily as others with use of applicable software.
407.3.2 Alphabetic Keys	Supports	Device has many keyboards and numeric keypad options in various configurations as well as access setting options with use of applicable software.



Criteria	Conformance Level	Remarks and Explanations
407.3.3 Numeric Keys	Supports	Device has many keyboards and numeric keypad options in various configurations as well as access setting options with use of applicable software.
407.4 Key Repeat	Supports	Device supports repeated button input. The delay before the key repeat feature is activated can be fixed at, or adjustable to, 2 seconds minimum with use of applicable software.
407.5 Timed Response	Does Not Support	The device hardware does not have applicable timed response.
407.6 Operation	Supports	Device supports at least one mode of operation operable with one hand and does not require tight grasping, pinching, or twisting of the wrist.
407.7 Tickets, Fare Cards, and Keycards	Not Applicable	Device is not a ticket, fare card, or keycard.
<b>407.8 Reach Height and Depth</b>	Heading cell – no response required	Heading cell – no response required
407.8.1 Vertical Reference Plane	Supports	This personal ICT device can be mounted on a vertical plane for side reach in a customizable position to support user requirements.
407.8.1.1 Vertical Plane for Side Reach	Supports	This personal ICT device can be mounted on a vertical plane for side reach in a customizable position to support user requirements.
407.8.1.2 Vertical Plane for Forward Reach	Supports	This personal ICT device can be mounted on a vertical plane for forward reach in a customizable position to support user requirements.
407.8.2 Side Reach	Supports	This personal ICT device can be mounted on a plane for side reach in a customizable position to support user requirements.

Criteria	Conformance Level	Remarks and Explanations
407.8.2.1 Unobstructed Side Reach	Supports	This personal ICT device can be mounted with unobstructed side reach in a customizable position to support user requirements.
407.8.2.2 Obstructed Side Reach	Supports	This personal ICT device can be mounted with obstructed side reach in a customizable position to support user requirements.
407.8.3 Forward Reach	Supports	This personal ICT device can be mounted with forward reach in a customizable position to support user requirements.
407.8.3.1 Unobstructed Forward Reach	Supports	This personal ICT device can be mounted with unobstructed forward reach in a customizable position to support user requirements.
407.8.3.2 Obstructed Forward Reach	Supports	This personal ICT device can be mounted with obstructed forward reach in a customizable position to support user requirements.
407.8.3.2.1 Operable Part Height for ICT with Obstructed Forward Reach	Supports	This portable personal ICT device can be mounted in any customizable position to support user access requirements.
407.8.3.2.2 Knee and Toe Space under ICT with Obstructed Forward Reach	Supports	This portable personal ICT device can be mounted in any customizable position to support user access requirements.
<b><u>408 Display Screens</u></b>	Heading cell – no response required	Heading cell – no response required
408.2 Visibility	Supports	This portable personal ICT device can be mounted in any customizable position to support user visibility of display screen.
408.3 Flashing	Supports	Device screen does not provide light flashes, with no more than three flashes in any one-second period.
<b><u>409 Status Indicators</u></b>	Heading cell – no response required	Heading cell – no response required

Criteria	Conformance Level	Remarks and Explanations
409.1 General	Supports	Device status indicators are discernible visually and by sound.
<b><u>410 Color Coding</u></b>	Heading cell – no response required	Heading cell – no response required
410.1 General	Supports	Device provides, color coding and also provides other means of conveying information, indicating an action, prompting a response, or distinguishing a visual element.
<b><u>411 Audible Signals</u></b>	Heading cell – no response required	Heading cell – no response required
411.1 General	Supports	Device provides audible signals or cues and also provides visual means of conveying information, indicating an action, or prompting a response.
<b><u>412 ICT with Two-Way Voice Communication</u></b>	Heading cell – no response required	Heading cell – no response required
<b>412.2 Volume Gain</b>	Heading cell – no response required	Heading cell – no response required
412.2.1 Volume Gain for Wireline Telephones	Not Applicable	Device is not a two-way voice communications piece of equipment.
412.2.2 Volume Gain for Non-Wireline ICT	Not Applicable	Device is not a two-way voice communications piece of equipment.
<b>412.3 Interference Reduction and Magnetic Coupling</b>	Heading cell – no response required	Heading cell – no response required
412.3.1 Wireless Handsets	Not Applicable	Device is not a two-way voice communications piece of equipment
412.3.2 Wireline Handsets	Not Applicable	Device is not a two-way voice communications piece of equipment
412.4 Digital Encoding of Speech	Not Applicable	Device is not a two-way voice communications piece of equipment
412.5 Real-Time Text Functionality	Reserved for future	Reserved for future
412.6 Caller ID	Not Applicable	Device is not a two-way voice communications piece of equipment
412.7 Video Communication	Not Applicable	Device is not a two-way voice communications piece of equipment but can support video communication while utilizing applicable software.

Criteria	Conformance Level	Remarks and Explanations
<b>412.8 Legacy TTY Support</b>	Heading cell – no response required	Heading cell – no response required
412.8.1 TTY Connectability	Not Applicable	Device is not a two-way voice communications piece of equipment
412.8.2 Voice and Hearing Carry Over	Not Applicable	Device is not a two-way voice communications piece of equipment
412.8.3 Signal Compatibility	Not Applicable	Device is not a two-way voice communications piece of equipment
412.8.4 Voice Mail and Other Messaging Systems	Not Applicable	Device is not a two-way voice communications piece of equipment
<b><u>413 Closed Caption Processing Technologies</u></b>	Heading cell – no response required	Heading cell – no response required
413.1.1 Decoding and Display of Closed Captions	Not Applicable	Device does not process video with synchronized audio or is a closed caption processing technology.
413.1.2 Pass-Through of Closed Caption Data	Not Applicable	Device does not process video with synchronized audio or is a closed caption processing technology.
<b><u>414 Audio Description Processing Technologies</u></b>	Heading cell – no response required	Heading cell – no response required
414.1.1 Digital Television Tuners	Not Applicable	Device does not digital television tuner or is an audio description processing technology.
414.1.2 Other ICT	Not Applicable	Device does not process video with synchronized audio or is an audio description processing technology.
<b><u>415 User Controls for Captions and Audio Descriptions</u></b>	Heading cell – no response required	Heading cell – no response required
415.1.1 Caption Controls	Not Applicable	Caption controls are not applicable to this device.
415.1.2 Audio Description Controls	Not Applicable	Audio Description controls are not applicable to this device.

## Chapter 5: Software

Notes:

Criteria	Conformance Level	Remarks and Explanations
501.1 Scope – Incorporation of WCAG 2.0 AA	See WCAG 2.x section	See information in WCAG 2.x section
<b><u>502 Interoperability with Assistive Technology</u></b>	Heading cell – no response required	Heading cell – no response required
502.2.1 User Control of Accessibility Features	Supports	Device supports User Control of Accessibility Features.
502.2.2 No Disruption of Accessibility Features	Supports	Device supports No Disruption of Accessibility Features.
<b>502.3 Accessibility Services</b>	Heading cell – no response required	Heading cell – no response required
502.3.1 Object Information	Not Applicable	Not Applicable
502.3.2 Modification of Object Information	Supports	Device States and properties that can be set by the user are capable of being set programmatically, including through assistive technology.
502.3.3 Row, Column, and Headers	Not Applicable	Device does not display data tables.
502.3.4 Values	Supports	Setting value(s), and any set or range of allowable values associated with an object are programmatically determinable.
502.3.5 Modification of Values	Supports	Values can be set by the user and are capable of being set programmatically, including through assistive technology.
502.3.6 Label Relationships	Supports	Relationships that a component has as a label for another component, or of being labeled by another component, are programmatically determinable.
502.3.7 Hierarchical Relationships	Supports	Hierarchical (parent-child) relationship that a component has as a container for, or being contained by, another component is programmatically determinable.
502.3.8 Text	Supports	The content of text objects, text attributes, and the boundary of text rendered to the screen, are programmatically determinable.
502.3.9 Modification of Text	Supports	Text that can be set by the user are

Criteria	Conformance Level	Remarks and Explanations
		capable of being set programmatically, including through assistive technology.
502.3.10 List of Actions	Supports	A list of all actions that can be executed on An object is programmatically determinable.
502.3.11 Actions on Objects	Supports	Device allows assistive technology to programmatically execute available actions on objects.
502.3.12 Focus Cursor	Supports	Device provides information and mechanisms necessary to track focus, text insertion point, and selection attributes of user interface components.
502.3.13 Modification of Focus Cursor	Supports	Text insertion point, and selection attributes that can be set by the user are capable of being set programmatically, including through the use of assistive technology.
502.3.14 Event Notification	Supports	Notification of events relevant to user interactions, including but not limited to, changes in the component's state(s), value, name, description, or boundary, are available to assistive technology.
502.4 Platform Accessibility Features	Supports	Device supports platform accessibility features.
<b><u>503 Applications</u></b>	Heading cell – no response required	Heading cell – no response required
503.2 User Preferences	Supports	Device permit user preferences from platform settings for color, contrast, font type, font size.
503.3 Alternative User Interfaces	Supports	Alternative User Interfaces.
<b><i>503.4 User Controls for Captions and Audio Description</i></b>	Heading cell – no response required	Heading cell – no response required
503.4.1 Caption Controls	Not Applicable	Caption controls are not applicable to this device.
503.4.2 Audio Description Controls	Not Applicable	Audio Description controls are not

Criteria	Conformance Level	Remarks and Explanations
		applicable to this device.
<b><u>504 Authoring Tools</u></b>	Heading cell – no response required	Heading cell – no response required
504.2 Content Creation or Editing (if not authoring tool, enter “not applicable”)	See <a href="#">WCAG 2.x</a> section	See information in WCAG 2.x section
504.2.1 Preservation of Information Provided for Accessibility in Format Conversion	Not Applicable	Device is not an authoring tool.
504.2.2 PDF Export	Not Applicable	Device is not an authoring tool.
504.3 Prompts	Not Applicable	Device is not an authoring tool.
504.4 Templates	Not Applicable	Device is not an authoring tool.

## Chapter 6: [Support Documentation and Services](#)

Notes:

Criteria	Conformance Level	Remarks and Explanations
<b><u>601.1 Scope</u></b>	Heading cell – no response required	Heading cell – no response required
<b><u>602 Support Documentation</u></b>	Heading cell – no response required	Heading cell – no response required
602.2 Accessibility and Compatibility Features	Supports	Device supports Accessibility and Compatibility Features.
602.3 Electronic Support Documentation	See <a href="#">WCAG 2.x</a> section	See information in WCAG 2.x section
602.4 Alternate Formats for Non-Electronic Support Documentation	Supports	Device supports Alternate Formats for Non-Electronic Support Documentation
<b><u>603 Support Services</u></b>	Heading cell – no response required	Heading cell – no response required
603.2 Information on Accessibility and Compatibility Features	Supports	Device supports Information on Accessibility and Compatibility Features.
603.3 Accommodation of Communication Needs	Supports	Device supports Accommodation of Communication Needs.

## PRC-Salttillo Legal Disclaimer

*This document is for informational purposes only. All product design and specification are subject to change.*

*Changes in the default settings of this application will have an impact on accessibility and may not adhere to accessibility best practices. In such cases, the responsibility for accessibility lies with the application owner.*

*The information provided in this report is true and correct to the best of our knowledge as of the date printed on the front page. The information provided for customer evaluation purposes only is subject without notice. The contents of this document do not constitute a legally binding claim. However, should you notice any issues, please contact PRC-Salttillo at [www.prc-salttillo.com/feedback](http://www.prc-salttillo.com/feedback) so that they may be evaluated and corrected, if necessary.*